

BEST[™] Compact SimulView[™] Ultrasound Imaging System

Compact System Features:

- Easily transportable—can be used with or without traditional stand used in conventional systems
- 21.5" touch screen technology
- Capable of supporting over 18 different probe configurations, the system is able to be used in almost every discipline



Best Medical International, Inc.
7643 Fullerton Road, Springfield, VA 22153 USA
tel 703 451 2378 800 336 4970
www.teambest.com www.nomos.com



A breakthrough in ultrasound technology—

BEST™ Sonalis® Imaging System

Featuring SimulView™ Ultrasound Technology



The Best™ Sonalis® Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView™ Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

Additional Features

- Patented probe design with simultaneous imaging of sagittal and transverse planes
- Longitudinal array provides for 140 mm length of view encompassing the bladder, prostate and perineum
- Superior HD Image resolution for improved implant accuracy, speed and physician confidence level
- Advanced drawing and editing tools which include user-defined line widths and colors for fiducial and anatomical markers
- Advanced modular software design provides for future upgrade path via in-house and independently developed technologies
- System converts from stand-based to desktop without losing any functionality
- Dual battery backup for an uninterrupted power supply in loss of power situations
- Multiple brachytherapy grids with user-definable grid symbols and colors
- Independent focal zones for sagittal and transverse images
- Independent frequency selection for sagittal and transverse Images
- New 21.5" touch screen technology

Advanced Keyboard Technology

Totally sealed, self-healing antibacterial keyboard with SensoFoil™ Technology



Sonalis[®] transducers

8L2A Linear Array

Applications: Arterial, Carotid, Vascular Access, Venous



12L5A Linear Array

Applications: Arterial, Breast, Carotid, Dialysis Access, Lung, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous



14L3 Linear Array

Applications: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous



15LW4 Linear Array

Applications: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous

Biopsy Kit Available



15LA Linear Array

Applications: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Testes, Thyroid, Vascular Access, Venous

Biopsy Kit Available



15L4A Linear Array

Applications: Arterial, Breast, Carotid, Dialysis Access, Lung, MSK, Neonatal Hip, Nerve Block, Ophthalmic, Thyroid, Vascular Access, Venous



16L5 Linear Array

Applications: Breast, Lung, MSK, Nerve Block, Vascular Access

VET Biopsy Kit Available



8V3 Phased Array

Applications: Cardiac



1 Best Drive Pittsburgh, PA 15202 USA

tel: 412 312 6700 800 70 NOMOS fax: 412 312 6701 www.nomos.com

We offer probes for the most common types of procedures and imaging—all with the enhanced image quality and high resolution that you would expect from the BEST name in ultrasounds!

4V2A Phased Array

Applications: Cardiac, FAST, TCD



5C2A Curved Array

Applications: Abdominal, FAST, Fetal Cardiac, MSK, OB/GYN, Renal, Thyroid, Visceral

Biopsy Kit Available



9MC3 Curved Array

Applications: Abdominal, Cardiac, Neonatal Head, Small Parts, Thyroid, Vascular Access



8EC4A Endocavity

Applications: OB/GYN, Prostate

Biopsy Kit Available



XY-BI-Plane Phased Array

Applications: Cardiac, Vascular, Lung



10EC4A Endocavity

Applications: OB/GYN, Prostate

Biopsy Kit Available



10BP4 Bi-Plane

Applications: Prostate



8BP4 Bi-Plane

Applications: Prostate



8TE3 Trans-esophageal

Applications: Motorized Adult Multiplane TEE Probe



Pedoff

Applications: Cardiac



16HL7 High Frequency Linear Array

Applications: MSK, Venous



Best™ **nomos**®
healthcare for everyone


TEAMBEST GLOBAL®



www.teambest.com

Best™ medical international

Best™ Total Solutions for Brachytherapy




Best™ Mick-Type Applicator



Best™ Kobold Henschke Tandem & Ovoid Applicator



Best™ Localization Needles with Best® I-125 Seeds





Best™ Esophageal Brachytherapy Applicator



Best™ Breast Double-Balloon Brachytherapy Applicator



Best™ Templates



Best™ Universal Chair/Table




Best™ medical international

Best™ Total Solutions for Brachytherapy

Best™ Fiducial Marker Kits



Best® Pd-103 Seeds



Best™ Fiducial Markers




Best Medical's Gold Fiducial Markers are available in a variety of diameters and lengths to meet your needs for soft tissue localization such as prostate, liver, pancreas, lung and breast.

- Best™ Gold Fiducial Markers are available in tissue-absorbable strands
- Standard strand consists of two markers at 2 cm apart center to center
- Special marker strands with custom spacing are also available
- Excellent visualization under X-ray, Fluoro, CT, MRI, Ultrasound and other imaging modalities
- Sold individually, stranded, and in pre-loaded, pre-sterilized needle sets
- Can be used for Brachytherapy, IMRT and IGRT

Best Medical is the only company making custom seeds and strands to your exact specifications—shipped within 24 hours, 7 days a week, sterile and non-sterile!

Best™ Brachytherapy Kit



5 Friction Cuffs

5 Implant Needles



5 Stylets with Hubs



5 Half Moon Buttons*

5 Single Leader Catheters**



5 Red Caps

*Half Moon Buttons are available in radiopaque (pictured) or clear nylon. **Catheter tubes are available in 4 colors (clear, blue, green or purple) with either radiopaque or clear nylon buttons.

Best Cure Global Healthcare Delivery System

15 Trillion U.S. Dollar Annual Global Economy



TeamBest Global companies (TBG) in collaboration with the Best Cure Foundation plan to establish 100s of Best Cure Global Healthcare Delivery of Proactive, Preventive, Primary Care and Multi-specialty Medical Centers to serve the needs of the global cancer population. They aim to promote Total Health™ (prevention, early detection, effective treatment and total cure) and practice Proactive Healthcare™ (full transparency on benefits, clinical outcome and costs) and reduce the cost of healthcare by 30 percent or more.

Best Cure Global Healthcare System is comprised of:

1. Best Cure Health System (non-profit)
2. Best Cure Insurance (non-profit)
3. Best Cure Clinical Research Institute (non-profit)
4. Best Medical Real Estate Investment Trust (for profit public company)
5. Best Medical Capital (for profit public company)
6. TBG companies (for profit public company)

All of the non-profit organizations will be under the Best Cure Foundation, a non-profit, global, private, non-governmental organization founded in 2007 and supported entirely by Krishnan Suthanthiran, President/Founder of TBG and BCF.

BestTM nomos[®]

1 Best Drive Pittsburgh, PA 15202 USA
Tel: 412 312 6700 800 70 NOMOS
Fax: 412 312 6701 www.nomos.com

BestTM medical international

7643 Fullerton Road Springfield, VA 22153 USA
Tel: 703 451 2378 800 336 4970
Fax: 703 451 5228 www.teambest.com

MSF-0019(R1)

The specifications contained in this brochure are subject to change.

TeamBest Global Companies ©2023